

Armedica Manufacturing Corp. has been in business for 10 years. Our founder and owner has been in the medical equipment business for over 40 years. He has been involved in many of the key design and manufacturing improvements in the treatment table industry. Our commitment to building solid products has helped us become one of the leading suppliers of treatment tables to the medical industry.

Made in the USA is not just a slogan at Armedica. We realize that designing, manufacturing and assembling all our products in our factory ensures that we have complete control over all the key processes. It also means that if you have special design requirements, we can make the changes and deliver a table that meets your specific needs.

We are an ISO13485 Certified Medical Device Manufacturer. In addition to being FDA compliant, we maintain a quality system that ensures the design and manufacture of every table meets global quality standards. That means a better and more reliable product for our customer.

We use the most reliable motors in the industry which are rated at 8000N (1798 lbs) of force. This costs us more to build each table, but it has proven to be the most trouble-free motor for the thousands of customers who own our tables.

Armedica Manufacturing Corp.
PO Box 880 | Greenwood, Arkansas 72936
800.701.5122 Toll Free | 479.996.3051 Fax | sales@armedicamfg.com

1 9 9 8 - 2 0 0 8
10
YEARS OF EXCELLENCE

MADE IN THE
USA



ISO13485
FM50089



Armedica
manufacturing corporation

The New **AM-BA SERIES** Bar Activated Hi-Lo Control

No more searching for the control switch!
This unique design makes the control easy for the clinician to use from anywhere around the table, yet out of the patient's way as they get on and off the table.



AM-BA300
Treatment Table
3-Section Top
With Elevating Center

**NEW Foot Activated
Hi-Lo Control!**

AM-BA150
Treatment Table
1-Section Top



AM-BA350
Treatment Table
3-Section Top
Non-elevating Center



*Also available in Elevating Center
Model: AM-BA300, shown on front cover*

AM-BA200
Treatment Table
2-Section Top



AM-BA400
Traction Table
4-Section Top



AM-BA234
Treatment Table
2-Section Top
Bariatric 34" Wide



*Also available in 40" BOBATH
Model: AM-BA240*

AM-BA450
Traction Table
4-Section Top
3-Piece Head Section



AM-BA250
Treatment Table
2-Section Top
3-Piece Head Section



AM-BA500
Treatment Table
3-Section Top
3-Piece Head Section
Elevating Center



AM-BA334
Treatment Table
34" Wide
3-Section Top
Non-elevating Center



Also available in 40" Model: AM-BA340

AM-BA550
Treatment Table
3-Section Top
3-Piece Head Section
Non-elevating Center



THE **NEW** BAR ACTIVATED HI-LO CONTROL AM-BA SERIES

Heavy duty steel tube Hi Lo control allows adjustment of the table height while keeping both hands on the patient and both feet on the floor. Control is designed to allow easy access from anywhere around the table without disrupting the treatment.

AM-BA Series Specifications Chart

- Four caster base with four inch (4") total lock casters standard on all AM-BA tables
- Optional Qualux® Ultra High Resiliency 2" Foam available for all AM-BA tables

ELEVATES (Blue)
LOWERS (Yellow)
BOTH (Green)
STATIONARY (White)

MODEL	LENGTH	WIDTH	CONFIGURATION	400lbs. CAP.	500lbs. CAP.	FACE HOLE
AM-BA150	76"	27"	[White] 1-SECTION TREATMENT TABLE	X		
AM-BA200	76"	27"	[Green] [White] 2-SECTION TREATMENT TABLE	X		X
AM-BA234	76"	34"	[Blue] [White] 2-SECTION BARIATRIC TREATMENT TABLE		X	X
AM-BA240	76"	40"	[Blue] [White] 2-SECTION BOBATH TREATMENT TABLE		X	
AM-BA250	76"	27"	[Yellow] [Green] [White] 4-SECTION TREATMENT TABLE	X		X
AM-BA300	76"	27"	[Green] [Blue] [Blue] 3-SECTION TREATMENT TABLE	X		X
AM-BA334	76"	34"	[Green] [White] [Blue] 3-SECTION BARIATRIC TREATMENT TABLE		X	X
AM-BA340	76"	40"	[Green] [Blue] [Blue] 3-SECTION BOBATH TREATMENT TABLE		X	X
AM-BA350	76"	27"	[Green] [White] [Blue] 3-SECTION TREATMENT TABLE	X		X
AM-BA400	76" *	27"	[Green] [White] [Blue] 4-SECTION TRACTION TABLE	X		X
AM-BA450	76" *	27"	[Yellow] [Green] [White] [Blue] 4-SECTION TRACTION TABLE	X		X
AM-BA500	76"	27"	[Green] [Blue] [Blue] 3-SECTION TREATMENT TABLE	X		X
AM-BA550	76"	27"	[Yellow] [White] [Blue] 3-SECTION TREATMENT TABLE	X		X

Available Colors

- Dove Gray
- Imperial Blue
- Tea Rose
- Cappuccino **
- Taupe
- Burgundy
- Blue Ridge
- Patina **
- Forest Green
- Black
- Merlot **

* Overall length with machine stand — 99 inches

** Decorator Colors have an ultra soft leather feel and soft hi-loft backing.