

Transfer & Position Aids

Medcare offers a full line of transfer and repositioning aids designed to assist in the mobility, handling, care, and transfer of patients or residents. Lifting, transferring and repositioning of patients and residents is the leading cause of back injury for caregivers.

For a needs assessment, free product demonstration or additional information, please contact us at 1-800-695-4479.

Motility Board

Ideal for assisted and non-assisted transfers, the Motility Board allows easy maneuvering from bed to chair, chair to wheelchair, etc. Lightweight, curved shape has a maximum safe working load of 600 lbs. Anti-slip pads provide maximum grip. An antimicrobial agent is built in for increased infection control.



Item Number:
255705

Dolphin Board

The Dolphin Toilet Transfer board is designed for assisted or independent transfer from wheelchair to toilet. The wheelchair cut-outs on the board allow transfer from either side of the toilet. Strong and lightweight, the board is rigid with no flex during transfer. Please note the board remains in position during use. The underside is treated with anti-slip material. It weighs 3.6 lbs. and has max user weight of 330 lbs.

Item Number:
255700



Pediturn

A soft low friction turning disc, with anti-slip material on outer surfaces. The Pediturn is positioned on the floor as an aid to turn a sufficient weight bearing person.


Item Numbers: 15.75" - 255301



Butterfly Board

Lightweight transfer board with unique removable nylon/polyurethane wing for controlled, low friction glide across the board. Suitable for all aspects of seated transfers wheelchair to chair, chair to bed, chair to car, etc. Suitable weight capacity of up to 265 lbs. The Butterfly Board has anti-slip pads on the underside for added security. Weight 3 lbs.

Item Numbers:
Board with medium wing - 255400
Board with large wing - 256989



Onward Turntable

Onward Turntable is a non-slip injection moulded turning disc. This turning disc is positioned on the floor as an aid to turn a sufficient weight bearing person. Disc is easy to disinfect between uses and is virtually maintenance free.


Item Number: 250006



Swivel Cushions

A low friction, circular swivel cushion is available in two standard sizes. The Swivel Cushion is made from comfortable quilted cotton on the top side and non-slip material on the underside. The swivel cushion can be positioned on a car seat, chair, or bed and gives assistance getting in or out of bed.


Item Numbers:
Small: 15.75" - 255300
Large: 19.75" - 255302



2 Move

The 2 Move range of products are designed for supine lateral transfers from bed to bed/stretcher/theatre table/x-ray/shower trolley. It is available in two widths and in adult and pediatric sizes.


Item Numbers:
15.75" x 33.5" - 256020,
nylon cover - 257001
19.75" x 33.5" - 256019,
nylon cover - 257002
15.75" x 67" - 254017,
nylon cover - 256995
19.75" x 67" - 254018,
nylon cover - 256996



Handling Belts

Available in Medium and Large sizes the Handling Belt enables better handling techniques while assisting transfer/repositioning. Adjustable and padded for increased patient comfort. An ideal aid when assisting with transfers using transfer boards or turning discs with one or two caregivers.

Item Numbers:
Premium Medium:
35" x 53" - 255304
Large: 45" x 65" - 255303
Standard Medium:
30" x 46" - 255306
Large: 40" x 56" - 255307



Premium Belt Shown

MedSheet and MedSheet Tube

The MedSheet is a repositioning device enables a caregiver to easily reposition patients/residents in a sitting or prone position. The MedSheet works to reduce friction by sliding on top of another MedSheet or folded over on top of itself. The MedSheet is available in sets of two sheets.

The MedSheet Tube is a loop of fabric with webbing handles used for repositioning a patient/resident on a bed.


Item Numbers:

MedSheet 67" x 34" – 400068

MedSheet 67" x 34" 10 pack – 400069

MedSheet Tube 40" x 32" – 400076

MedSheet Tube 76" x 33" – 400095



MedSheet



MedSheet Tube



MedSheets can be used in the following situations:

- Moving a fallen patient/resident from a spot not accessible with a lift.
- Repositioning a prone patient/resident that has slid down the bed.
- Repositioning a prone patient/resident to a seated position.

Features

- Designed with resident and caregiver safety in mind.
- Silicone coating provides a low-friction surface which allows the caregiver to easily slide the resident without lifting.
- Durable reinforced “box x” stitching provides safe handles for sliding.
- Rip-stop technology built in.



While the MedSheet is a versatile tool, it should not be used as a substitute for vertical lifting equipment such as a Medicare Lift or Stand. Use of the MedSheet for lifting patients could result in serious injury to a caregiver and/or patient.

One-Way-Glide Anti-Slip Sheets

Used to help maintain a good sitting position in the wheelchair/day chair/bed. Available in 3 sizes, the One-Way-Glide has been designed to move in the direction of the arrow only, and then helps prevent the patient from slipping forward. Machine wash at 170° F. Hang to dry or tumble on low 104° F.

Item Number:

Small chair size: 17" x 14.5" – 255200

Chair size: 17.75" x 17.75" – 255201

Bed size: 31.5" x 23.5" – 255202




At Medcare Products, our mission is to provide the highest-quality patient lifting, transferring and positioning equipment, along with the highest standard of training and service in the industry.

Antimicrobial Protection


Toileting, incontinence, wound care and simple everyday use make slings and belts very susceptible to microbe growth. It is also common practice for slings and belts to be used multiple times on multiple people during the day before they are washed. Medicare Slings and Belts are the first and *only* to be manufactured with Microban antimicrobial product protection.

Microban® protection provides continuous antimicrobial cleaning that inhibits the growth of microbes such as bacteria, mold and mildew that can cause stains, odors and product deterioration. Built-in during the manufacturing process, Microban antimicrobial protection will not wash or wear off, providing protection for the useful lifetime of the product and keeping the slings and belts cleaner between cleanings. Test results prove its effectiveness against bacteria.




Bacteria Performance

These samples were tested for efficacy against bacteria using the TM 147 test method. Both samples were inoculated with bacteria and then photographed after 24 hours. The top sample was manufactured with Microban protection and the bottom sample was untreated. The sample with built-in Microban protection showed a significant inhibition zone surrounding the treated sling material. The untreated material showed no inhibition zone.



Medicare Sling Material
Untreated



Medicare Sling Material
Treated with Microban Protection



Reposition Aids

Transfer & Reposition Aids

Accessories

One-Way-Glide Anti-Slip Sheets

Used to help maintain a good sitting position in the wheelchair/day chair/bed. Available in 3 sizes, the One-Way-Glide has been designed to move in the direction of the arrow only, and then helps prevent the patient from slipping forward. Machine wash at 170° F. Hang to dry or tumble on low 104° F.

Item Number:

Small chair size: 17" x 14.5" – 255200

Chair size: 17.75" x 17.75" – 255201


Bed size: 31.5" x 23.5" – 255202



Multiglide Gloves

The Multiglide glove is an effective aid in reducing the friction when handling patients. The glove is carefully glided under the client and is an effective tool for repositioning or for checking therapeutic seating surfaces.

Item Number: 250201



MedSheet

The MedSheet is composed of two nylon sheets which allow for easy movement of a patient during repositioning.

Item Number: 400068



MedSheet Tube

The MedSheet Tube is a loop of fabric with webbing handles used for repositioning a patient in a bed.

Item Number: 400095



Body Support Wedge

These are manufactured from flame-retardant foam and are supplied complete with a removable, waterproof, brushed polyurethane cover, which is also flame-retardant. The wedge provides a positive non-slip action when used against cotton bed linen. The wedge also gives a firm support for the patient.

Item Numbers:

Wedge with cover – 256984

Additional cover – 256979



Washable Cotton Overpads

Two sizes are available. Protective cotton over pads, for use on the One-Way-Glides and 4 Way Glide. Also for general use to protect day chairs, wheelchairs and beds from the problems associated with incontinence.

Item Numbers:

17" x 24" – 255203

34" x 36" –

255204



At Medcare Products, our mission is to provide the highest-quality patient lifting, transferring and positioning equipment, along with the highest standard of training and service in the industry.

www.medcarelifts.com | 800-695-4479

151 East Cliff Road, Burnsville, MN 55337

©2008 Medcare, Inc. Specifications are subject to change without notice.

All products referred herein are registered trademarks or trademarks of their respective holders.